

## REHAB SHOWER COMMODES

Traditionally, rehab shower commodes have been constructed of steel or stainless steel. These two materials had their drawbacks. Steel rusts, stainless steel is expensive and both are heavy. Ezee Life™ commodes are constructed of high quality aluminum that won't rust, is lighter than steel and lower in cost than stainless steel.


**TWO CASTER  
OPTIONS:  
STANDARD OR LONG  
STEMMED CASTERS  
STANDARD STEM PART # SSC  
LONG STEM PART # LSC**



New locking flip up arms.



Optional 4 Way Seat



Anti tippers  
on all models




Height adjustable  
on all models



|  |  |  |
|--|--|--|
| <p>Open Front Seat &amp; Insert</p> <p>Standard on all sizes</p> | <p>Closed Front Seat &amp; Top</p> <p>Available on 180, 185 &amp; 186 Models</p> | <p>4-Way Open Seat</p> <p>Available on 180</p> |
|--|--|--|

Removable snap on seat (all models)



### STANDARD FEATURES:

- Lightweight aluminum **NO RUST**
- Footrest, adjustable height & removable
- Removable insert cushion for showering or toileting
- Quick Release axles (self propelled models)
- Four 5" locking casters
- Slide on pail, removable
- Flip back arms
- Maintenance free tires
- Calf strap included
- Molded handrims
- Super soft seat
- Easier rolling
- Seat belt
- All adjustments are tool-less
- 1 Year full warranty
- 2 Year warranty on frame

| COMMODOE MODELS                                      |            |
|--|------------|
| Model Codes  | Seat Width |
| 150  | 15"        |
| 180  | 18"        |
| 185  | 20"        |
| 186  | 22"        |
| Model numbers indicate the seat size of the commode. |            |

| OPTIONAL EQUIPMENT AND FEATURES                             |                       |                         |
|---|-----------------------|-------------------------|
| Option Codes  | Option Name           | Availability            |
| 4W  | 4 Way seat            | All 18" Commodes        |
| CSF18   | 18" Closed Front Seat | All 18" Commodes        |
| CSF 20  | 20" Closed Front Seat | All 20" Commodes        |
| CSF 22  | 22" Closed Front Seat | All 22" Commodes        |
| 22  | 22" Rear Wheels       | 15" Commodes            |
| 24  | 24" Rear Wheels       | 18", 20" & 22" Commodes |
| *Add the option codes above to the model numbers to select  |                       |                         |
| *Open Front Seat and 5" Rear Casters Standard on All Models |                       |                         |

| SPECIFICATIONS   | 150, 150-22        | 180, 180-4W, 180-24, 180-4W-24 | 185 Wide, 185-24   | 186 X-Wide, 186-24  |
|--|--------------------|--------------------------------|--------------------|---------------------|
| Frame Material   | Aluminum           | Aluminum                       | Aluminum           | Aluminum            |
| Weight Capacity  | 286 lbs (130 kgs)  | 286 lbs (130 kgs)              | 385 lbs (175 kgs)  | 440 lbs (200 kgs)   |
| Rear Wheels  | 5" locking casters | 5" locking casters             | 5" locking casters | 5" locking casters  |
| Optional Rear Wheels                                       | 22" (150-22)       | 24" (180-24)                   | 24" (185-24)       | 24" (186-24)        |
| Seat Width   | 15"                | 18"                            | 20"                | 22"                 |
| Seat Depth   | 15"                | 17"                            | 17"                | 17"                 |
| Foot Rest Hanger Angle                                     | 60°                | 60°                            | 60°                | 60°                 |
| Foot Rest Length   | 14" to 20"         | 14" to 20"                     | 14" to 20"         | 14" to 20"          |
| Back Height  | 18"                | 18"                            | 18"                | 18"                 |
| Product Weight   | 26.4 lbs (12 kgs)  | 28.4 lbs (12.9 kgs)            | 29 lbs (13.2 kgs)  | 32.3 lbs (14.7 kgs) |
| Arm Height   | 11" (from seat)    | 11" (from seat)                | 11" (from seat)    | 11" (from seat)     |
| Arm Length   | 13"                | 13"                            | 13"                | 13"                 |
| Arm Pad Length   | 9"                 | 9"                             | 9"                 | 9"                  |
| Tilt Range Posterior                                       | n/a                | n/a                            | n/a                | n/a                 |
| Rear Internal Width  | 12"                | 15"                            | 17"                | 19"                 |
| Seat Heights (Standard casters)                            | 19.5" or 20.5"     | 19.5" or 20.5"                 | 19.5" or 20.5"     | 19.5" or 20.5"      |
| Seat Heights (Long stem casters)                           | 21.5" or 22.5"     | 21.5" or 22.5"                 | 21.5" or 22.5"     | 21.5" or 22.5"      |
| Seat Clearance Height <sup>2</sup> (Standard caster stems) | 15.5" or 16.5"     | 15.5" or 16.5"                 | 15.5" or 16.5"     | 15.5" or 16.5"      |
| Seat Clearance Height <sup>2</sup> (Long caster stems)     | 17.5" or 18.5"     | 17.5" or 18.5"                 | 17.5" or 18.5"     | 17.5" or 18.5"      |

| DIMENSIONS WITH 5" REAR CASTERS | 150 | 180 | 185 Wide | 186 X-Wide |
|---------------------------------|-----|-----|----------|------------|
| Overall Width                   | 19" | 22" | 24"      | 26"        |
| Overall Length <sup>1</sup>     | 37" | 37" | 37"      | 37"        |

| DIMENSIONS WITH 24" REAR WHEELS | 150 | 180 | 185 Wide | 186 X-Wide |
|---------------------------------|-----|-----|----------|------------|
| Overall Width                   | 23" | 26" | 28"      | 30"        |
| Overall Length <sup>1</sup>     | 42" | 42" | 42"      | 42"        |

1) Length with all parts extended (wheels highest height, anti-tip bars fully extended, footrests in place)

2) Seat clearance is the measurement from the bottom of the pail rail to the floor.