

How it works

Open your vehicle's door and fold the **XL-SEAT** down. Once transferred onto the **XL-SEAT**, use the retractable support arm for safety and balance while it takes you up to the driver's seat.

Features

- Keep your original seat
- All seat adjustment functions remain operational
- Seat back reclines completely
- Removable transfer board to allow easy access for others


Specifications

- Maximum capacity: 275 lbs
- Transfer board size: 15¾" X 17"
- Lowers to the height of the rocker panel (floor of the vehicle) and can be stopped at any height, according to your needs
- Seat covered with vinyl surface and slight cushioning for easy transfer

Installation

The **XL-SEAT** is intended for the driver and front passenger positions in Minivans, Pickup Trucks, SUV's and some Crossovers with no structural modification to your vehicle.

The way the **XL-SEAT** was designed makes it the first choice when looking for a mobility solution that offers simplicity, comfort and space. It is preferred by many because of its compact style and design.



Optional
HI-LIFT



XL-SEAT

Motorised transfer board

Options

We offer two solutions to get your wheelchair into your minivan via the sliding door: the **HI-LIFT** or the **SPEEDY-LIFT**.

HI-LIFT - Optional wheelchair lift

The **HI-LIFT** is specially designed as an option, to work as a pair with the **XL-SEAT**.

Features

- Fully motorised lifting and lowering functions
- Motorised in and out swing mechanism
- Compact design takes up very little space

Specifications

- Maximum capacity: 75lbs
- Designed for manual foldable or rigid frame wheelchairs

Due to the diversity of mobility solutions available to you, please contact your local mobility dealer for a personalized evaluation and overview of options tailored to your needs.

Our products qualify as adaptive equipment and are eligible for reimbursement under the terms of Mobility programs offered by various vehicle manufacturers.

Safety

We take our customers' safety seriously. We are dedicated to manufacturing quality mobility equipment. Our products undergo a strict certification process. Certification of our products involves different tests such as: crash testing, cycle testing and weight capacity testing. We are proud to offer a line of products that exceed the industry's standards.

Dealer :

ADAPT SOLUTIONS

QUALITY MOBILITY EQUIPMENT

145 DAMASE-BRETON | SAINT-LAMBERT, QC, G0S 2W0

P 866 641 0419 | F 418 889 9838 | O 418 889 0419

www.adaptsolutions.ca



creation : www.noirdstudio.com



www.adaptsolutions.ca