

## Go-Go Scooter Comparison Chart

Scooter	Traveller		Traveller Plus		Ultra X	
<b>Configuration:</b>	3 Wheeled	4 Wheeled	3 Wheeled	4 Wheeled	3 Wheeled	4 Wheeled
<b>Drive Type:</b>	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive
<b>Maximum Speed:</b>	Up to 6 kph (4 mph)	Up to 6 kph (4 mph)	Up to 7.2 kph (4.5 mph)	Up to 7.2 kph (4.5 mph)	Up to 6 kph (4 mph)	Up to 6 kph (4 mph)
<b>Max Capacity:</b>	136 kg (300 lb)	136 kg (300 lb)	136 kg (300 lb)	136 kg (300 lb)	118 kg (260 lb)	118 kg (260 lb)
<b>Distance on Full Battery*:</b>	Up to 11.5 km (7.2 mi)	Up to 11.5 km (7.2 mi)	Up to 19.8 km (12.4 mi)	Up to 21.3 km (13.3 mi)	Up to 12.8 km (8 miles)	Up to 12.8 km (8 miles)
<b>Battery:</b>	Sealed 12 A/H battery (optional 18 A/H avail.)	Sealed 12 A/H battery (optional 18 A/H avail.)	2 x 18 A/H Sealed	2 x 18 A/H Sealed	2 x 12 A/H	2 x 12 A/H
<b>Turning Radius:</b>	84 cm (33")	112 cm (44")	86.4 (34")	116.2 cm (45.75")	83.8 cm (33")	112 cm (44")
<b>Weight:</b>	34 kg (75 lb)	35.7 kg (78.5 lb)	38.2 kg (84 lb)	39.3 kg (86.5 lb)	32.7 kg (72 lb)	35.7 kg (78.5 lb)
<b>Ground Clearance:</b>	5 cm (2")	5 cm (2")	7.6 cm (3")	7.6 cm (3")	5 cm (2")	5 cm (2")
<b>Suspension:</b>	None	None	None	None	None	None

\* Scooter range dependent on user weight, battery health, tire condition, and local terrain.

Scooter	Sport		LX		ES 2	Folding
<b>Configuration:</b>	3 Wheeled	4 Wheeled	3 Wheeled	4 Wheeled	3 Wheeled	4 Wheeled (Dual front)
<b>Drive Type:</b>	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive
<b>Maximum Speed:</b>	Up to 6 kph (4 mph)	Up to 6 kph (4 mph)	Up to 7 kph (4.4 mph)	Up to 7 kph (4.4 mph)	Up to 7 kph (4.4 mph)	Up to 6 kph (4 mph)
<b>Max Capacity:</b>	147.7 kg (325 lb)	147.7 kg (325 lb)	136 kg (300 lb)	136 kg (300 lb)	113.6 cm (250 lb)	113.6 cm (250 lb)
<b>Distance on Full Battery*:</b>	Up to 20.3 km (12.7 mi)	Up to 19.2 km (12 miles)	Up to 15.1 km (9.44 mi)	Up to 16.8 km (10.5 mi)	Up to 9.6 km (6 mi)	Up to 14.9 km (9.3 mi)
<b>Battery:</b>	2 x 18 A/H Sealed	2 x 18 A/H Sealed	Sealed 12 A/H battery (optional 18 A/H avail.)	Sealed 12 A/H battery (optional 18 A/H avail.)	2 x 10 A/H Sealed	Sealed 12 A/H battery (optional 11.5 AH lithium-ion)
<b>Turning Radius:</b>	91.4 cm (36")	132 cm (52")	84 cm (33")	113.7 cm (44.75")	81.3 cm (32")	81.3 cm (32")
<b>Weight:</b>	37.8 kg (83.2 lb)	45.1 kg (99.2 lb)	38.8 kg (85.4 lb)	39.7 kg (87.4 lb)	33.6 kg (74 lb)	31.5 kg (69.4 lb)
<b>Ground Clearance:</b>	7.6 cm (3")	7.6 cm (3")	3.8 cm (1.5")	3.8 cm (1.5")	5 cm (2")	5.8 cm (2.3")
<b>Suspension:</b>	None	None	Comfort-Trac independent suspension	Comfort-Trac independent suspension	None	None

\* Scooter range dependent on user weight, battery health, tire condition, and local terrain.