

Pride Scooter Comparison Chart

Scooter	Victory 10		Maxima		Victory Sport
Configuration:	3 Wheeled	4 Wheeled	3 Wheeled	4 Wheeled	4 Wheeled
Drive Type:	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive
Maximum Speed*:	Up to 8.4 kph (5.2 mph)	Up to 8.4 kph (5.2 mph)	Up to 8.4 kph (5.2 mph)	Up to 8.4 kph (5.2 mph)	Up to 12.8 kph (8 mph)
Max Capacity:	181 kg (400 lb)	181 kg (400 lb)	227.2 kg (500 lb)	227.2 kg (500 lb)	159 kg (350 lb)
Distance on Full Battery*:	Up to 24.8 km (15.5 mi)	Up to 24.8 km (15.5 mi)	Up to 24.6 km (15.4 mi)	Up to 24.6 km (15.4 mi)	Up to 17.6 km (11 miles)
Battery:	(2) U1 or (2) 40AH	(2) U1 or (2) 40AH	(2) NF-22	(2) NF-22	(2) 40AH
Turning Radius:	106.7 cm (42")	137.1 cm (54")	125.2 cm (49.3")	137.1 cm (54")	131.4 cm (51.75")
Weight:	59.8 kg (131.5 lb)	60.2 kg (132.5 lb)	84.5 kg (186 lb)	87.3 kg (192 lb)	68.2 kg (150 lb)
Ground Clearance:	7.6 cm (3")	7.6 cm (3")	8.2 cm (3.25")	8.2 cm (3.25")	5.7 cm (2.25")
Suspension:	None	None	Rear	Front & Rear	Full Suspension

* Scooter speed and range dependent on user weight, battery health, tire condition, and local terrain.

Scooter	Pursuit	Pursuit XL	Pursuit Sport (36v)	Sport Rider	
Configuration:	4 Wheeled	4 Wheeled	4 Wheeled	3 Wheel / Single Seat	3 Wheel / Dual Seat
Drive Type:	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive	Rear Wheel Drive
Maximum Speed*:	Up to 11.6 kph (7.25 mph)	Up to 13.6 kph (8.5 mph)	14.4 kph (9 mph)	Up to 15 kph (9.4 mph)	Up to 12 kph (7.5 mph)
Max Capacity:	181 kg (400 lb)	181 kg (400 lb)	181 kg (400 lb)	181 kg (400 lb)	181 kg (400 lb)
Distance on Full Battery*:	Up to 22.4 km (14 mi)	G27: Up to 47.2 km G24: Up to 35.4 km	Up to 36 km (22.5 mi)	Up to 50.4 km (31.25 mi)	Up to 52.5 km (32.8 mi)
Battery:	(2) NF-22	(2) Grp. 27 or (2) Grp. 24	(3) 12-volt, deep cycle, 55 AH	(2) 12-volt, deep cycle (75 AH or 100 AH)	(2) 12-volt, deep cycle (75 AH or 100 AH)
Turning Radius:	170.8 cm (67.25")	209.5 cm (82.5")	210.8 cm (83")	158 cm (62.2")	158 cm (62.2")
Weight:	99.3 kg (218.5 lb)	118.1 kg (259.9 lb)	123.2 kg (271 lb)	115.2 kg (253.5 lb)	116.9 kg (257.2 lb)
Ground Clearance:	8.9 cm (3.5")	13.3 cm (5.25") (midframe)	6.3 cm (2.5")	12.7 cm (5")	12.7 cm (5")
Suspension:	Full Suspension	Full Suspension	Full Suspension	Full Suspension	Full Suspension

* Scooter speed and range dependent on user weight, battery health, tire condition, and local terrain.